



## SWITCH 2.0 HPC

### THE BENCHMARK

Designed for daily use in theatre, intensive care and emergency vehicles, the reusable ERKA. Infusion Cuffs offer sophisticated details and efficient handling. The material used – a polyester fabric and nylon – is ideally suited for efficient pressure buildup followed by uniform emptying of the infusion bag. The infusion cuff has a quality manometer with an ergonomically shaped ball – for rapid pumping. The constant influx of infusion solutions, even in emergency situations, is therefore guaranteed at all times. The ERKA infusion cuffs SWITCH 2.0 HPC are available in two sizes and suitable for all standard infusion bags with a volume of 500 ml or 1,000 ml. Suitable for the MR environment – on a par with all the other innovative products belonging to the ERKA SWITCH 2.0 brand. Transparent – thanks to the transparent net, the bag and its content can be clearly seen. Quality – hook and high quality manometer are included.

#### ARE THERE DIFFERENT SIZES AVAILABLE?

The ERKA infusion cuffs SWITCH 2.0 HPC are available in two sizes: 500 and 1,000 ml.

#### FOR WHICH INFUSION BAGS ARE THE ERKA HPC SUITABLE?

ERKA high pressure cuffs are suitable for all standard infusion bags.

#### ARE ERKA HPC MR-CONDITIONAL?

As part of the innovative ERKA. Switch 2.0 product family, the pressure infusion cuff Switch 2.0 HPC is specially adapted to the particular conditions of the MRI environment.

#### LONG-LASTING

The mesh fabric of the cuff is durable and resilient. The transparent mesh material allows for a clear view of the bag and its contents.

#### HYGIENIC

The removable outer sleeve of the cuff is washable and easy to disinfect.

#### SKIN-FRIENDLY

All materials used are dermatologically safe and latex-free.

Switch 2.0 HPC

500 ml

Size

-

Reference

822.00000.

Accessories

Size

Reference

Spare Parts

Pump Ball grey	000.20900
Pump Ball with Reducing Valve, grey	000.21100
Measuring Vessel + Hook	000.30201
Cover 500 ml	822.30000
Cover 1.000 ml	832.30000
Bladder 500 ml	821.30000
Bladder 1.000 ml	831.30000

## Legal Notes and Data Protection

ERKA.  
Kallmeyer Medizintechnik GmbH & Co. KG  
Im Farchet 15  
83646 Bad Tölz

Geschäftsführer: Björn Kallmeyer

HRA 75980, Amtsgericht München  
Umsatzsteuer-ID: DE 812 833 893

Kontakt:  
Telefon: 0 80 41 / 80 09-0  
Telefax: 0 80 41 / 80 09-39  
eMail: [made@erka.org](mailto:made@erka.org)  
<http://www.erka.org>

Verantwortlich:  
Björn Kallmeyer

## Signatur / Stempel

SKIN-FRIENDLY

